



## BLOOD COLLECTION AND SUBMISSION FORM

### Owner/Breeder Information

---

**Name:** \_\_\_\_\_ **Date:** \_\_\_\_\_

---

**Address:** \_\_\_\_\_ **City:** \_\_\_\_\_ **State:** \_\_\_\_\_ **Zip Code:** \_\_\_\_\_

---

**Phone:** (     )                      **Cell :** (     )                      **Email:** \_\_\_\_\_

---

**Name of Owner/Business to be printed on certificate:** \_\_\_\_\_

---

**Comments:** \_\_\_\_\_

---

**Need more kits?**    **Yes**    **No**                      **How many?** \_\_\_\_\_

---


	Bird's Name or ID	Bird Species
1		
2		
3		
4		
5		

### Instructions For Collecting Blood Sample

1. Wash your hands thoroughly before and after handling each bird. (Clean clippers with alcohol before and after each bird)
2. Clearly label each blood card with the bird's name and/or ID number and the species of the bird. Use the same Bird ID and/or name on this form so that we can match them to the correct bird later.
3. Once you have safely restrained your bird, cut the tip of the birds nail using clean clippers. Rub the nail to help the blood flow. Use a cotton ball to clean off the first drop of blood. Dab 2-3 drops of blood on the blood card. (Use styptic powder or flour to stop the bleeding if necessary)
4. Allow the blood card to air dry completely. Place the blood card in a plastic bag.
5. Complete this form and put it in an envelope with the blood samples and payment and send it to:

**JP Bird Sexing**  
 5403 West 22<sup>nd</sup> Court  
 Hialeah, FL 33016

**Pricing**  
 1-4 samples \$15.00 each  
 5 or more samples \$13.00 each



**Choose a Payment Option:**

\_\_\_ Paypal Transaction ID: \_\_\_\_\_

\_\_\_ Check or Money Order made out to JP Bird Sexing Amount Enclosed: \_\_\_\_\_